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Note intended for the approved laboratories

Datum/Date: 05/08/2013
DIENSTMITTEILUNG / OFFICE CIRCULAR
Nr. / nr. KO/LABO/1097702

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Subject : Detection of *Cronobacter sakazakii* in dried milk

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The method laid down in the UE regulation 2073/2005/EC for the detection of *Cronobacter sakazakii* in dried milk is the ISO/TS 22964 method.

Yet the prescribed isolation medium (ESIA-agar) and some confirmation tests do not enable to differentiate *Cronobacter sakazakii* from *Cronobacter malonaticus*. In this context, the National Reference Laboratory for Food Microbiology has verified in a study which commercial isolation media and confirmation tests do enable to differentiate those two strains.

Besides the ESIA agar, a second isolation medium can be chosen between :

- ChromAgar *E. sakazakii* (Biotrading) – powder: both strains show blue colonies. The *Cronobacter sakazakii* colonies are less intense than the *Cronobacter malonaticus* colonies.
- Rapid Sakazakii agar (Biorad) – ready to use: *Cronobacter sakazakii* show green-grey colonies. *Cronobacter malonaticus*

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The following tests are recommended as confirmation tests (on visual reading only):

- Rapid ID32 E: on **visual** reading, differences can be observed with following tests: MNT (sodium malonate test), PLE (palatinose), ONAG (β -D-glucosaminide). The MNT and ONAG tests are positive for *Cronobacter malonaticus* and negative for *Cronobacter sakazakii*, and the opposite for the PLE test.
- API 20E: on **visual** reading, differences can be observed with following tests: VP (Voges Proskauer), INO (Inositol). These tests are negative for *Cronobacter malonaticus* and positive for *Cronobacter sakazakii*

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Crystal ID does not differentiate between *Cronobacter malonaticus* and *Cronobacter sakazakii*.

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We ask the laboratories to implement these recommendations in their method.

Generaldirektor

Director general

ir. Geert De Poorter (unterz./sgd)